





DISPOSABLE STERILE OPHTHALMIC SURGICAL KNIFE






Slit Knife

No.	Model	Bend angle $\alpha 1 \pm 5^\circ$	Edge width $W \pm 0.05$ (mm)	Blade thickness t (mm)
1	CJY-01-1.0	45°	1.0	0.18-0.22
2	CJY-01-1.1		1.1	
3	CJY-01-1.2		1.2	
4	CJY-01-1.3		1.3	
5	CJY-01-1.4		1.4	
6	CJY-01-1.6		1.6	
7	CJY-01-1.8		1.8	
8	CJY-01-2.0		2.0	
9	CJY-01-2.1		2.1	
10	CJY-01-2.2		2.2	
11	CJY-01-2.3		2.3	
12	CJY-01-2.4		2.4	
13	CJY-01-2.5		2.5	
14	CJY-01-2.6		2.6	
15	CJY-01-2.7		2.7	
16	CJY-01-2.8		2.8	
17	CJY-01-2.9		2.9	
18	CJY-01-3.0		3.0	
19	CJY-01-3.2		3.2	
20	CJY-01-3.5		3.5	
20	CJY-01-1.0 (Z)	Straight	1.0	
22	CJY-01-2.8 (Z)	Straight	2.8	
23	CJY-01-3.0 (Z)	Straight	3.0	

Crescent Knife




No.	Model	Bend angle $\alpha 1 \pm 5^\circ$	Edge width $W \pm 0.05$ (mm)	Blade thickness t (mm)
1	CJY-02-2.0	Angled 60° , bevel-up	2.0	0.18-0.22
2	CJY-02-2.3	Angled 60° , bevel-up	2.3	
3	CJY-02- <u>2.3</u>	Angled 60° , bevel-down	<u>2.3</u>	
4	CJY-02-2.0 (Z)	Straight	2.0	
5	CJY-02-2.3 (Z)	Straight	2.3	

R Slit Knife



No.	Model	Bend angle $\alpha 1 \pm 5^\circ$	Edge width $W \pm 0.05$ (mm)	Blade thickness t (mm)
1	CJY-03-2.6	45°	2.6	0.18-0.22
2	CJY-03-3.2		3.2	
3	CJY-03-3.5		3.5	
4	CJY-03-3.7		3.7	
5	CJY-03-4.1		4.1	
6	CJY-03-5.2		5.2	
7	CJY-03-5.5		5.5	




MVR Knife

No.	Model	Bend angle $\alpha \pm 5^\circ$	Edge width $W \pm 0.05$ (mm)	Pipe diameter $\phi \pm 0.02$ (mm)
1	CJY-04-19g	45°	1.5	1.0
2	CJY-04-20g		1.2	0.88
3	CJY-04-23g		0.9	0.53
4	CJY-04-19g (Z)	Straight	1.5	1.0
5	CJY-04-20g (Z)	Straight	1.2	0.88
6	CJY-04-23g (Z)	Straight	0.9	0.53



Straight Knife

No.	Model	Blade angle α	Edge width $W \pm 0.05$ (mm)	Blade thickness t (mm)
1	CJY-05-0°	0°	1.6	0.16-0.22
2	CJY-05-15°	15°		
3	CJY-05-22.5°	22.5°		
4	CJY-05-30°	30°		
5	CJY-05-45°	45°		



Pterygium knife

No.	Model	Bend angle α	Edge width $W \pm 0.05$ (mm)	Blade thickness t (mm)
1	CJY-06-2.4	Straight	2.4	0.20-0.32

Product Function

Designed for ophthalmic microsurgery

Product Features

Exclusive automation, quality guarantee, high-level sharpness, above industry standard

The high-level sharpness ensures less strength in slitting and moving in operation so that the incision edge is smooth and the eyeball moves slightly.

[Production License No.]
Su Food and Drug Manufacturer Production No. 20160082


[Number of Registration Certificate]
CFDA Applicable 20162041338

About Us

Founded in 2003, Jiajie Medical is located in Changzhou City, Jiangsu Province, China. Our company is devoted to R&D, manufacturing, marketing and service of disposable medical instruments. Our products include disposable surgical knives for ophthalmology and other disposable medical knives (annular knife, linear knife, overmold knife, endoscopic knife, cambered knife and so on). We have won the trust and praise from customers by the stable quality.

In the philosophy of "Quality First", "Integrity Based" and "Constant Improvement", our company is determined to become a leader in the industry. We actively join the international market and accept the advanced information and technology from the industry so that we can offer the most professional service to customers.

Our domestic customers includes more than fifty well-known enterprises such as Weigao Group, Panther Healthcare, Reach Surgical, Frankenman, Suzhou Touchstone, Brightness, Jiangxi Langhe, Kangdi, Waston, etc.



No.8-21, Shunyuan Road, Xinbei, Changzhou, Jiangsu, China.

Cutter knife sales Tel / +86-519-8822 0928

Ophthalmic knife sales Tel / +86-519-8518 6899

Fax / +86-519-8822 0968

E-mail / ivy@jsjjajie.com

